



## Hygienic Orthodontic Acrylic



One of the finest on the market today. Can be used for repair or making appliances. Cure time 6 to 8 minutes.

**1 lb. Standard Package**  
1 lb. Powder and 8 oz. of Liquid.

**5 lbs. Lab Package**  
5 lbs. Powder, 1 quart of Liquid and 1—8 oz. of Liquid.



Clear	Pink	Description	Each
0131-000	0131-001	1 lb. Standard Package	<b>\$56.90</b>
0131-003	0131-004	5 lb. Lab Package	<b>\$229.00</b>
0131-A	0131-A	1 lb. Powder Only	<b>\$41.90</b>
0131-C	0131-X	8 oz. Liquid Only	<b>\$20.90</b>
0131-M	0131-P	1 Quart Liquid Only	<b>\$41.90</b>

## MASEL® Ultracryl II Orthodontic Acrylic



**New Improved Formula!**

For superior strength and hardness for fog-free, transparent appliances.

Acrylic strength is essential when making orthodontic removable appliances. **Ultra-cryl**, unlike conventional orthodontic retainer acrylic, can withstand the constant occlusal stress demanded of functional appliances and TMJ appliances. Not only is it highly crack-resistant, but it has also been formulated to set slowly and have little flow, providing ample time for precision detailing of appliances. This quick-cure, cold polymerizing acrylic mixes well, as the powder quickly absorbs the liquid, avoiding slumping and running - so often found with other acrylics. This lets you fabricate the appliance with a minimum bulk and maximum strength. Available in light pink and clear shades, both of which are resistant against oral acidity. Best results-with pressure pot.

- Powder Only—1 lb. or 5 lb. Can
- Liquid Only—1 Quart (32 oz)



Clear	Pink	Description	Each
0131-LA2	0131-LA2	1 lb. Powder Only	<b>\$41.90</b>
0131-HA2	0131-KA2	1 qt. Liquid Only	<b>\$39.90</b>
0131-NA2	0131-NA2	5 lb. Powder Only	<b>\$138.00</b>

## MASEL® Electrik™ Color Concentrate



Ideal conclusion for patients preferring colored accessories during treatment. Mix 10 drops of color concentrate into clear monomer. Add more drops to achieve deeper color. 2 oz. (60cc).

Radiant Red	2032-103
Bold Blue	2032-104
Electrik Yellow	2032-105
Black Blaze	2032-106
Emerald Green	2032-107
Vivid Violet	2032-108

*Color concentrate to be added to clear ultracryl liquid or other methyl-methacrylate base clear acrylic*

**\$18.90 ea.**

## Electrik™ Neon Ultracryl



Fun and flashy retainers make patient compliance skyrocket! Mix and match to create distinctive designs. Polymer and monomer combine to make a strong, durable retainer. Includes directions. 2 oz. (50 g).

Neon Orange	2032-462
Neon Red	2032-463
Neon Glow	2032-464

*Use neon polymer with clear monomer. Then mix and apply like regular acrylic.*

**\$18.90 ea.**

## Electrik™ Glitter



Anything that can make kids wear their retainers is a plus. Add it to your powder or sprinkle into acrylic before setting. Squeeze bottle with fine nozzle eliminates waste, allows you to create individual design. 2 oz. (60cc).

Silver	2032-300
Gold	2032-311
Red	2032-312
Blue	2032-313
Confetti	2032-320

*Mix 5 squeezes of glitter with 2 oz. of regular polymer. Then use like regular acrylic.*

**\$14.90 ea.**

## Super Sep™ Separating Medium

Superior alginate based separating solution seals pores on model surface, making acrylic appliances remove easily and safely. Ready-to-use. 1 liter.

# 0131-Z  
**\$22.95 ea.**

